


# Colorectal Anastomosis Leak Testers

**Expedite and standardize the 'bubble test' for colorectal anastomoses**



**ANASTOMOSIS LEAK TESTER - LP**  
inflates to 30 cm water above  
intraoperative pressure when abdomen  
is insufflated to 15 mm Hg

White for LAPAROSCOPIC surgery  
**Catalog No. ALT-50**

## Anastomosis Leak Testers:

- Completely occlude the anus
- Simplify rectal inflation
- Signal when recommended pressure\* is reached
- Prevent over-distension
- Color coded
- Sterile, compact and disposable

Patent Pending

\*Wheeler JMD, Gilbert JM. Controlled intraoperative water testing of left-sided colorectal anastomoses: are ileostomies avoidable? *Ann R Coll Surg Engl* 1999;81:105-8

Black for OPEN surgery  
**Catalog No. ALT-30**



**ANASTOMOSIS LEAK TESTER - OP**  
inflates to 30 cm water

**TfS** Tools for Surgery

**We make it simple, by design.**

Tools for Surgery, LLC  
1339 Stony Brook Road ■ Stony Brook, New York 11790  
PH (631) 444-4448 ■ FX (631) 444-5330

[www.toolsforsurgery.com](http://www.toolsforsurgery.com)