

# SIREN SGT™

36 French orogastric tube for bariatric and gastric surgery

## SIREN SGT™

Catalog No. SGT-36

Supplied clean, not sterile  
5 units per box

Pinch clamp  
simplifies opening  
and closing

Integral suction  
regulating valve  
limits suction pressure  
to a safe and  
effective level

Prominent numerical  
markings indicate  
depth of insertion

Blunt tip, smooth  
surface, and flexibility  
facilitate safe insertion

Multiple side holes  
facilitate gastric  
decompression



U.S. Patent No. 11,006,957

**TfS** Tools  
for  
Surgery

[www.toolsforsurgery.com](http://www.toolsforsurgery.com)

Tools for Surgery, LLC  
8 Technology Drive, #100 • East Setauket, NY 11733  
P: 631-444-4448 • F: 631-689-2943  
info@toolsforsurgery.com