

Sohn's Dilators®

for Standardized Anal Dilatation for anal fissure

AUTOCLAVABLE

Catalog No. SOD-03
Set of three dilators



- Standardized Anal Dilatation is an effective treatment for chronic anal fissure, with fewer complications and a lower rate of incontinence compared to lateral internal sphincterotomy.
- Success rates for anal dilatation to 40 mm for 6 minutes are 83% to 88%.
- Sohn's Dilators® are simpler and more economical than pneumatic balloons.

Sohn et al. Precise anorectal dilatation - its role in the therapy of anal fissures. Dis Colon Rectum 1992;35:322-327

Sohn N, Weinstein MA. Anal dilatation for anal fissures. Seminars in Colon & Rectal Surgery 1997;8(1):17-23

Pinsk et al. The long-term effect of Standardized Anal Dilatation for chronic anal fissure on anal continence. Ann Coloproctol 2021;37(2):115-119

Abe et al. Long-term efficacy and safety of controlled manual anal dilatation in the treatment of chronic anal fissures: a single-center observational study. J Anus Rectum Colon 2023;7(4):250-257



www.toolsforsurgery.com/SOD

Tools for Surgery, LLC
8 Technology Drive, #100 • East Setauket, NY 11733
P: 631-444-4448 • F: 631-689-2943
info@toolsforsurgery.com